

RIVER[®]

RIVER[®]

Coronary angioplasty balloon catheter

NEW, LOW PROFILE ANGIOPLASTY BALLOON CATHETER:

- *Dedicated for narrow and tortuous vessels*
- *Extraordinary flexibility*
- *Excellent pushability and trackability*

Caution:

This publication could include technical or other inaccuracies or typographical errors. Changes are periodically added to the information herein; these changes will be incorporated in the new editions of the publication.
Balton Ltd. may make improvements and/or changes in the product(s) described in this publication at any time.

01/2013



ul. Modlińska 294
03-152 Warszawa, Poland

www.balton.pl
balton@balton.pl

tel. +48 (22) 597 44 00
fax +48 (22) 597 44 44



RIVER[®]

Coronary angioplasty balloon catheter



Technical data :

Balloon material	nylon
Nominal pressure	8 atm
Catheter coating	hydrophilic
Tip entry profile	0.017" (0.43 mm)
Crossing profile	0.025" (0.64 mm)
Proximal shaft diameter	1.8F (0.60 mm)
Distal shaft diameter	2.5F (0.84 mm)
Pushability and trackability	excellent
Inflation and deflation time	very short
Compatible guidewire	max. 0.014"
Compatible guiding catheter OD/ID	5F / 0.056"
Catheter length	140 cm
Delivery system	Rapid Exchange

Reference: Balloon diameter / length

CBR1.25x10140	1.25 / 10 mm
CBR1.25x15140	1.25 / 15 mm
CBR1.25x20140	1.25 / 20 mm
CBR1.25x25140	1.25 / 25 mm
CBR1.25x30140	1.25 / 30 mm
CBR1.25x35140	1.25 / 35 mm
CBR1.25x40140	1.25 / 40 mm
CBR1.50x10140	1.50 / 10 mm
CBR1.50x15140	1.50 / 15 mm
CBR1.50x20140	1.50 / 20 mm
CBR1.50x25140	1.50 / 25 mm
CBR1.50x30140	1.50 / 30 mm
CBR1.50x35140	1.50 / 35 mm
CBR1.50x40140	1.50 / 40 mm
CBR1.75x10140	1.75 / 10 mm
CBR1.75x15140	1.75 / 15 mm
CBR1.75x20140	1.75 / 20 mm
CBR1.75x25140	1.75 / 25 mm
CBR1.75x30140	1.75 / 30 mm
CBR1.75x35140	1.75 / 35 mm
CBR1.75x40140	1.75 / 40 mm
CBR2.00x10140	2.00 / 10 mm
CBR2.00x15140	2.00 / 15 mm
CBR2.00x20140	2.00 / 20 mm
CBR2.00x25140	2.00 / 25 mm
CBR2.00x30140	2.00 / 30 mm
CBR2.00x35140	2.00 / 35 mm
CBR2.00x40140	2.00 / 40 mm
CBR2.25x10140	2.25 / 10 mm
CBR2.25x15140	2.25 / 15 mm
CBR2.25x20140	2.25 / 20 mm
CBR2.25x25140	2.25 / 25 mm
CBR2.25x30140	2.25 / 30 mm
CBR2.25x35140	2.25 / 35 mm
CBR2.25x40140	2.25 / 40 mm
CBR2.50x10140	2.50 / 10 mm
CBR2.50x15140	2.50 / 15 mm
CBR2.50x20140	2.50 / 20 mm
CBR2.50x25140	2.50 / 25 mm
CBR2.50x30140	2.50 / 30 mm
CBR2.50x35140	2.50 / 35 mm
CBR2.50x40140	2.50 / 40 mm

Reference: Balloon diameter / length

CBR2.75x10140	2.75 / 10 mm
CBR2.75x15140	2.75 / 15 mm
CBR2.75x20140	2.75 / 20 mm
CBR2.75x25140	2.75 / 25 mm
CBR2.75x30140	2.75 / 30 mm
CBR2.75x35140	2.75 / 35 mm
CBR2.75x40140	2.75 / 40 mm
CBR3.00x10140	3.00 / 10 mm
CBR3.00x15140	3.00 / 15 mm
CBR3.00x20140	3.00 / 20 mm
CBR3.00x25140	3.00 / 25 mm
CBR3.00x30140	3.00 / 30 mm
CBR3.00x35140	3.00 / 35 mm
CBR3.00x40140	3.00 / 40 mm
CBR3.25x10140	3.25 / 10 mm
CBR3.25x15140	3.25 / 15 mm
CBR3.25x20140	3.25 / 20 mm
CBR3.25x25140	3.25 / 25 mm
CBR3.25x30140	3.25 / 30 mm
CBR3.25x35140	3.25 / 35 mm
CBR3.25x40140	3.25 / 40 mm
CBR3.50x10140	3.50 / 10 mm
CBR3.50x15140	3.50 / 15 mm
CBR3.50x20140	3.50 / 20 mm
CBR3.50x25140	3.50 / 25 mm
CBR3.50x30140	3.50 / 30 mm
CBR3.50x35140	3.50 / 35 mm
CBR3.50x40140	3.50 / 40 mm
CBR3.75x10140	3.75 / 10 mm
CBR3.75x15140	3.75 / 15 mm
CBR3.75x20140	3.75 / 20 mm
CBR3.75x25140	3.75 / 25 mm
CBR3.75x30140	3.75 / 30 mm
CBR3.75x35140	3.75 / 35 mm
CBR3.75x40140	3.75 / 40 mm
CBR4.00x10140	4.00 / 10 mm
CBR4.00x15140	4.00 / 15 mm
CBR4.00x20140	4.00 / 20 mm
CBR4.00x25140	4.00 / 25 mm
CBR4.00x30140	4.00 / 30 mm
CBR4.00x35140	4.00 / 35 mm
CBR4.00x40140	4.00 / 40 mm